

## Enzyme Immunoassay for Lipoxin A<sub>4</sub>

For Research Use Only

---

### INTRODUCTION

---

Lipoxin A<sub>4</sub> (LXA<sub>4</sub>) is a biologically active lipoxygenase interaction product derived from arachidonic acid. Arachidonic acid is first oxygenated by 15-lipoxygenase to form 15-HETE which is converted by 5-lipoxygenase and epoxide hydrase to generate LXA<sub>4</sub>. LXA<sub>4</sub> stimulates leukocyte chemotaxis without aggregation and inhibits natural killer cell cytotoxicity. It also provokes contraction of parenchymal strips and stimulates microvascular changes. Recent findings indicate that LXA<sub>4</sub> inhibits leukocyte-dependent inflammation. Determination of LXA<sub>4</sub> level may provide new understanding of the role of LXA<sub>4</sub> in basic cellular reactions and in pathophysiology of inflammation and other disease processes.

---

### PRINCIPLES OF PROCEDURE

---

This is an ELISA (Enzyme-Linked ImmunoSorbent Assay) for the quantitative analysis of Lipoxin A<sub>4</sub> levels in biological fluid. This test kit operates on the basis of competition between the enzyme conjugate and the LXA<sub>4</sub> in the sample for a limited number of antibody binding sites.

The sample or standard solution is first added to the microplate. Next, the diluted enzyme conjugate is added and the mixture is shaken and incubated at room temperature for one hour. During the incubation, competition for binding sites is taking place. The plate is then washed removing all the unbound material. The bound enzyme conjugate is detected by the addition of substrate that generates an optimal color after 30 minutes. Quantitative test results may be obtained by measuring and comparing the absorbance reading of the wells of the samples against the standards with a microplate reader at 650 nm. The extent of color development is inversely proportional to the amount of LXA<sub>4</sub> in the sample or standard. For example, the absence of LXA<sub>4</sub> in the sample will result in a bright blue color, whereas the presence of LXA<sub>4</sub> will result in decreased or no color development.

---

### MATERIALS PROVIDED

---

Component	Description	Volume	Storage
<b>EIA Buffer</b>	Buffer used to dilute the Conjugate and LXA <sub>4</sub> Standards.	30 mL	4°C
<b>10x Wash Buffer</b>	Buffer used to wash the plate prior to color development.	20 mL	4°C
<b>TMB Substrate</b>	TMB substrate used for color development.	20 mL	4°C
<b>5x Extraction Buffer</b>	Buffer used to dilute extracted and non-extracted samples.	30 mL	4°C
<b>LXA<sub>4</sub>-HRP Conjugate</b>	Lyophilized LXA <sub>4</sub> horseradish peroxidase conjugate.	2 vials	-20°C
<b>LXA<sub>4</sub> Standard</b>	1 µg/mL LXA <sub>4</sub> standard solution.	100 µL	4°C
<b>Coated Plate</b>	96-well microplate coated with rat anti-LXA <sub>4</sub> antibody.	1 plate	4°C

---

### MATERIALS NEEDED BUT NOT PROVIDED

---

1. Microplate reader with a 450 nm or 650 nm filter
2. Adjustable micropipettes (10 – 1000 µL) and tips
3. Deionized water

4. Plate cover or plastic film
5. Test tubes
6. 1 N HCl (optional)

---

**EXTRACTION MATERIALS**

---

1. 1 N HCl
2. Methanol
3. C<sub>18</sub> Sep-Pak® Light Column (Waters® Corporation #23501)
4. Hexane
5. Methyl Formate
6. Nitrogen Gas
7. Centrifuge

---

**STORAGE**

---

1. Store the components of this kit at the temperatures specified on the labels.
2. Unopened reagents are stable until the indicated kit expiration date.
3. Concentrated, reconstituted conjugate has a shelf life of at least two weeks when stored properly.
4. Desiccant bag must remain in foil pouch with unused strips. Keep pouch sealed when not in use to maintain a dry environment. Remove excess air before sealing.

---

**WARNINGS AND PRECAUTIONS**

---

1. Use aseptic technique when opening and dispensing reagents.
2. This kit is designed to work properly as provided and instructed. Additions, deletions or substitutions to the procedure or reagents are not recommended, as they may be detrimental to the assay.

---

**PROCEDURAL NOTES**

---

1. Ensure that the conjugate is completely reconstituted before use. Each vial, when reconstituted, provides sufficient reagent to perform 64 wells (8 strips). If more than 64 assays are to be run, reconstitute both vials and pool the reconstituted conjugate.
2. To minimize errors in absorbance measurements due to handling, wipe the exterior bottom of the microplate wells with a lint-free paper towel prior to inserting into the plate reader.

---

**SAMPLE PREPARATION**

---

Plasma, serum, urine and tissue culture need to be extracted.

---

**EXTRACTION PROTOCOL**

---

1. Dilute 100  $\mu$ L of sample with 200  $\mu$ L of methanol, and then dilute that volume with 1.5 mL of water.
2. Acidify 1 mL of biological fluid to pH 3.5 with 1 N HCl (about 150  $\mu$ L for plasma) and vortex.
3. For tissue, homogenize it in ethanol (5 mL/g). Centrifuge the homogenate for five minutes. Collect the supernatant in a clean tube. Dilute 1 mL of the supernatant with 5 mL of deionized water and acidify to pH 3.5 with 1 N HCl.
4. Precondition the C<sub>18</sub> Sep-Pak® Light column by washing the column with 2 mL of methanol followed by 2 mL of deionized water.
5. Apply the above sample into the column and adjust the flow rate to 1 mL per minute. Reducing the flow rate to 0.5 mL per minute may increase extraction efficiencies.
6. Wash the column with 5 mL of water followed by 5 mL of hexane.
7. The LXA<sub>4</sub> is eluted by 2 mL of methyl formate.

8. Evaporate the methyl formate eluate with a stream of nitrogen gas.
9. Resuspend the residue with 1 mL of diluted Extraction Buffer. The residue may be dissolved in less than 1 mL if the concentration is suspected to be low (<0.1 ng per mL).

---

## REAGENT PREPARATION

---

1. **5x Extraction Buffer:** Dilute the appropriate amount to 1x with deionized water prior to use.
  2. **10x Wash Buffer:** Add 20 mL of 10x Wash Buffer to 180 mL of deionized water prior to use.
  3. **LXA<sub>4</sub>-HRP Conjugate:** Reconstitute by adding 75  $\mu$ L of deionized water to each vial. Dilute 110  $\mu$ L of Conjugate into 5.5 mL total volume of EIA Buffer.
- 

## STANDARD CURVE PREPARATION

---

The LXA<sub>4</sub> Standard is provided as a 1  $\mu$ g/mL stock solution. Use the following tables to dilute a set of standard stock solutions and construct an eight-point standard curve.

**Table 1:** Standard Stock Preparation

Standard	LXA <sub>4</sub> Conc. (ng/mL)	Vol. of EIA Buffer ( $\mu$ L)	Transfer Vol. ( $\mu$ L)	Final Vol. ( $\mu$ L)
A	1000	-	Provided	80
B	20	980	20 $\mu$ L of A	800
C	2	1800	200 $\mu$ L of B	1800
D	0.2	1800	200 $\mu$ L of C	2000

**Table 2:** Standard Curve Preparation

Standard	LXA <sub>4</sub> Conc. (ng/mL)	Vol. of EIA Buffer ( $\mu$ L)	Vol. of Stock C ( $\mu$ L)	Vol. of Stock D ( $\mu$ L)
S <sub>0</sub>	0	1000	-	-
S <sub>1</sub>	0.02	900	-	100
S <sub>2</sub>	0.04	800	-	200
S <sub>3</sub>	0.1	500	-	500
S <sub>4</sub>	0.2	-	-	1000
S <sub>5</sub>	0.4	800	200	-
S <sub>6</sub>	0.8	600	400	-
S <sub>7</sub>	2.0	-	1000	-

---

## ASSAY PROCEDURE

---

1. Add 50  $\mu$ L of Standards or Samples (may require diluting) to the corresponding wells on the microplate in duplicate. See **Scheme I** for a sample plate layout.
2. Add 50  $\mu$ L of diluted LXA<sub>4</sub>-HRP Conjugate to each well. Incubate at room temperature for one hour.
3. Wash the plate three times with 300  $\mu$ L of diluted Wash Buffer per well. Wash 5 times if using an automated plate washer.
4. Add 150  $\mu$ L of TMB Substrate to each well. Incubate at room temperature for 30 minutes.
5. Read the plate at 650 nm.

Alternately, the color reaction can be stopped after 10-15 minutes by adding 50  $\mu$ L of 1 N HCl and read at 450 nm.


**NOTE:** If accounting for substrate background, use 2 wells as blanks (BLK) with only 150  $\mu$ L TMB Substrate in the wells. Subtract the average of these absorbance values from the absorbance values of the wells being assayed.

**Scheme I: Sample Plate Layout**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	S0	S0	U1	U1	U9	U9	U17	U17	U25	U25	U33	U33
<b>B</b>	S1	S1	U2	U2	U10	U10	U18	U18	U26	U26	U34	U34
<b>C</b>	S2	S2	U3	U3	U11	U11	U19	U19	U27	U27	U35	U35
<b>D</b>	S3	S3	U4	U4	U12	U12	U20	U20	U28	U28	U36	U36
<b>E</b>	S4	S4	U5	U5	U13	U13	U21	U21	U29	U29	U37	U37
<b>F</b>	S5	S5	U6	U6	U14	U14	U22	U22	U30	U30	U38	U38
<b>G</b>	S6	S6	U7	U7	U15	U15	U23	U23	U31	U31	U39	U39
<b>H</b>	S7	S7	U8	U8	U16	U16	U24	U24	U32	U32	BLK	BLK

**CALCULATIONS**

1. After the substrate background has been subtracted from all absorbance values, average all of your duplicate well absorbance values.
2. The average of your two S<sub>0</sub> values is now your B<sub>0</sub> value. (S<sub>1</sub> now becomes B<sub>1</sub>, etc.)
3. Next, find the percent of maximal binding (%B/B<sub>0</sub> value). To do this, divide the averages of each standard absorbance value (now known as B<sub>1</sub> through B<sub>7</sub>) by the B<sub>0</sub> absorbance value and multiply by 100 to achieve percentages.
4. Graph your standard curve by plotting the %B/B<sub>0</sub> for each standard concentration on the y-axis against concentration on the x-axis. Draw a curve by using a curve-fitting routine (i.e. 4-parameter or linear regression).
5. Divide the averages of each sample absorbance value by the B<sub>0</sub> value and multiply by 100 to achieve percentages.
6. Using the standard curve, the concentration of each sample can be determined by comparing the %B/B<sub>0</sub> of each sample to the corresponding concentration of standard.
7. If the samples were diluted, the concentration determined from the standard curve must be multiplied by the dilution factor.

**Figure 1: Typical Standard Curve**

---

**CROSS REACTIVITY**

---

Lipoxin A <sub>4</sub>	100.0%	12-HETE	<0.10%
15-epi- Lipoxin A <sub>4</sub>	24%	Leukotriene B <sub>4</sub>	<0.01%
5(S)6(R)-Di-HETE	5.00%	Leukotriene C <sub>4</sub>	<0.01%
Lipoxin B <sub>4</sub>	1.00%	Leukotriene D <sub>4</sub>	<0.01%
15-HETE	0.10%	Leukotriene E <sub>4</sub>	<0.01%
5-HETE	<0.10%		

---

**DISCLAIMER**

---

This information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Oxford Biomedical Research, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See catalog for additional terms and conditions of sale.

---

**ORDERING INFORMATION**

---

For additional kits or a complete catalog please call 800-692-4633.

---

**TECHNICAL SUPPORT**

---

If you need technical information or assistance with assay procedures, call our Technical Support Department at 800-692-4633 or 248-852-8815. Our staff will be happy to answer your questions about this or any other product in the Oxford Biomedical line.

---

**GUARANTEE AND LIMITATION OF REMEDY**

---

Oxford Biomedical Research, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the materials in this kit, except that these materials and this kit will meet our specifications at the time of delivery. Buyer's remedy and Oxford Biomedical Research, Inc.'s sole liability hereunder is limited to, at Oxford Biomedical Research, Inc.'s option, refund of the purchase price of, or the replacement of, material that does not meet our specification. By acceptance of our products, Buyer indemnifies and holds Oxford Biomedical Research, Inc. harmless against, assumes all liability for the consequence of its use or misuse by the Buyer, its employees, or others. Said refund or replacement is conditioned on Buyer notifying Oxford Biomedical Research, Inc. within thirty (30) days of the receipt of product. Failure of Buyer to give said notice within thirty (30) days of receipt of product shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

Oxford Biomedical Research, Inc.  
P.O. Box 522  
Oxford, MI 48371 U.S.A.

Orders: 800-692-4633  
Technical Service: 248-852-8815  
Fax: 248-852-4466  
E-mail: [info@oxfordbiomed.com](mailto:info@oxfordbiomed.com)

Made in the U.S.A.