

## Enzyme Immunoassay for Human Myeloperoxidase

For Research Use Only

---

### INTRODUCTION

---

Human myeloperoxidase (MPO) is a hemoprotein composed of two heme-containing heavy subunits of 53 kDa and of two light subunits of 15 kDa. MPO is stored in the granules of neutrophils and monocytes, and is released in response to leukocyte activation. MPO acts as a catalyst in the production of hypochlorous acid (HOCl), a powerful oxidant produced from chloride ion (Cl<sup>-</sup>) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). Although MPO is important to the immune system, it also plays a role in several inflammatory conditions.

Elevated MPO levels have been linked to coronary artery diseases. High concentrations of MPO and its oxidation products are found in areas of plaque deposits present in atherosclerosis. MPO also affects cholesterol levels through the oxidation of LDL and HDL.

As a specific marker for polymorphonuclear cells (PMNs), MPO is released extracellularly via degranulation after PMN activation. Thus, extracellular MPO is an index of PMN activation in inflammatory processes that lead to disease pathology, such as cardiovascular disease. Elevated levels of MPO present in serum and plasma indicate a greater risk of myocardial injury and unstable plaque formations. Therefore, MPO levels can be used to predict potentially negative cardiac events.

Antibodies to MPO are present in human plasma, and have been characterized as part of the anti-neutrophil cytoplasmic antibodies (ANCA) set of autoantibodies.

---

### PRINCIPLES OF PROCEDURE

---

The Human MPO EIA Kit is a sandwich ELISA (Enzyme-Linked Immunosorbent Assay) that measures MPO in plasma/serum samples. The MPO in the sample interacts with and binds to the anti-MPO antibody coated on the wells. The plate is then washed to remove all unbound materials and anti-MPO, conjugated with HRP (horseradish peroxidase), is added to the plate wells. The HRP conjugate binds to the MPO that is bound to the anti-MPO of the plate. After washing the plate again, TMB is added, which will produce a blue hue. A darker blue color corresponds to a higher concentration of MPO. After reaction has occurred with the conjugate a stop solution is added, producing a yellow hue. The absorbance is measured at 450 nm.

---

### MATERIALS PROVIDED

---

Component	Description	Volume	Storage	Cat. No.
<b>5x EIA Buffer</b>	Buffer used to dilute the conjugate, standard, and samples.	50 mL	4°C	MP14a
<b>20x Wash Buffer</b>	Buffer used to wash the plate.	50 mL	4°C	MP14b
<b>Stop Solution</b>	0.3 M sulfuric acid for stopping the reaction.	12 mL	4°C	MP14c
<b>TMB</b>	TMB substrate used for color development.	12 mL	4°C	MP14d
<b>Anti-MPO-HRP</b>	Anti-MPO horseradish peroxidase conjugate.	200 µL	4°C	MP14e
<b>MPO Standard</b>	Lyophilized human myeloperoxidase.	1 vial	4°C	MP14f
<b>Coated Plate</b>	96-well microplate coated with an anti-Hu MPO antibody.	1 plate	4°C	MP14g

---

**MATERIALS NEEDED BUT NOT PROVIDED**

---

1. Microplate reader with a 450 nm filter
2. Adjustable micropipettes (10 – 1000  $\mu$ L) and tips
3. Deionized water
4. Vortex
5. Plate cover or film
6. Plate shaker/rotator

---

**STORAGE**

---

1. Store the components of this kit at the temperatures specified on the labels.
2. Unopened reagents are stable until the indicated kit expiration date.
3. Store diluted reagents at 4°C and use within one week.
4. Desiccant packet must remain in the pouch with unused strips. Keep pouch sealed when not in use to maintain a dry environment. Remove excess air before sealing.

---

**WARNINGS AND PRECAUTIONS**

---

1. Use aseptic technique when opening and dispensing reagents.
2. The MPO Standard has been purified from human source material. Necessary precautions should be taken when handling, treating it as potentially infectious.
3. This kit is designed to work properly as provided and instructed. Additions, deletions or substitutions to the procedure or reagents are not recommended, as they may be detrimental to the assay.

---

**PROCEDURAL NOTES**

---

1. Allow kit components to come to room temperature before use.
2. The HRP conjugate is most stable in its concentrated form. Dilute only the volume necessary for the amount of strips currently being used.
3. After reconstitution of the MPO Standard, aliquots should be made and frozen at -80°C. Avoid multiple freeze/thaw cycles.
4. To minimize errors in absorbance measurements due to handling, wipe the exterior bottom of the microplate wells with a lint-free paper towel prior to inserting into the plate reader.

---

**SAMPLE COLLECTION AND STORAGE**

---

This ELISA is meant for the measurement of myeloperoxidase in human serum and plasma.

1. Draw blood into a serum separator tube for serum or a heparin tube for plasma and mix gently.
2. Centrifuge whole blood within 6 hours of draw at 3000 x g for 10 minutes at 4°C.
3. Remove the serum supernatant. Samples can be stored at 4°C for 24 hours. For longer storage, samples should be stored at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

---

**SAMPLE PREPARATION**

---

Thaw samples completely. Dilute the samples in 1x EIA Buffer. 1:4 and 1:10 are recommended.

---

**REAGENT PREPARATION**

---

1. **5x EIA Buffer:** Dilute the appropriate amount to 1x with deionized water prior to use.
2. **20x Wash Buffer:** Add 50 mL of 20x Wash Buffer to 950 mL of deionized water prior to use.
3. **Anti-MPO-HRP:** Dilute 1:100 - 100  $\mu$ L of Conjugate into 9.9 mL total volume of 1x EIA Buffer.
4. **MPO Standard:** Reconstitute with 1 mL deionized water.

---

## STANDARD CURVE PREPARATION

---

The MPO Standard is provided as a lyophilized powder, and needs to be reconstituted with 1 mL of deionized water prior to use. Use the following table to construct a six-point standard curve.

**Table 1:** Standard Curve Preparation

Standard	MPO Conc. (ng/mL)	Vol. of 1x EIA Buffer ( $\mu$ L)	Transfer Volume ( $\mu$ L)	Transfer Source	Final Volume ( $\mu$ L)
S5	30	538	62	Stock	300
S4	15	300	300	S5	300
S3	7.5	300	300	S4	300
S2	3.75	300	300	S3	300
S1	1.875	300	300	S2	600
B0	0	300	-	-	300

---

## ASSAY PROCEDURE

---

1. Add 100  $\mu$ L of Standards and Samples to the corresponding wells on the microplate in duplicate. See **Scheme I** for a sample plate layout. Cover and incubate at room temperature for one hour.
2. Wash the plate four times with 300  $\mu$ L of 1x Wash Buffer per well.
3. Add 100  $\mu$ L of diluted Anti-MPO-HRP to each well. Cover and incubate at room temperature for 20 minutes.
4. Wash the plate four times with 300  $\mu$ L of 1x Wash Buffer per well.
5. Add 100  $\mu$ L of TMB to each well. Place covered plate on a plate shaker and incubate at room temperature for 10 minutes.
6. Add 100  $\mu$ L of Stop Solution to each well.
7. Read the plate at 450 nm.

**Scheme I:** Sample Plate Layout

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	B0	B0	U3	U3	U11	U11	U19	U19	U27	U27	U35	U35
<b>B</b>	S1	S1	U4	U4	U12	U12	U20	U20	U28	U28	U36	U36
<b>C</b>	S2	S2	U5	U5	U13	U13	U21	U21	U29	U29	U37	U37
<b>D</b>	S3	S3	U6	U6	U14	U14	U22	U22	U30	U30	U38	U38
<b>E</b>	S4	S4	U7	U7	U15	U15	U23	U23	U31	U31	U39	U39
<b>F</b>	S5	S5	U8	U8	U16	U16	U24	U24	U32	U32	U40	U40
<b>G</b>	U1	U1	U9	U9	U17	U17	U25	U25	U33	U33	U41	U41
<b>H</b>	U2	U2	U10	U10	U18	U18	U26	U26	U34	U34	U42	U42

---

## CALCULATIONS


---

1. Subtract the average background absorption (B0) from each of the mean absorbance values ( $A_{450}$ ) for the standards and samples.
2. Graph the standard curve by plotting the  $A_{450}$  for each standard concentration on the y-axis against the MPO concentration on the x-axis. Draw the curve by using a linear regression curve-fitting routine. Most programs can do this for you. See **Figure 1** for a sample standard curve.
3. Using the standard curve and the resulting linear regression equation, the concentration of MPO in each sample can be determined using the  $A_{450}$  of each sample.

Linear Regression Equation:  $y = mx + b$   
 Where:  $y = A_{450}$      $x = \text{MPO Concentration}$   
 $m = \text{slope}$      $b = \text{y-axis intercept}$

- If the samples were diluted, the concentration determined from the standard curve must be multiplied by the dilution factor.

**Figure 1: Typical Standard Curve**




**NOTES**

- The normal ranges for serum samples and heparin–plasma samples are 2.5–100 ng/mL and 37-124 ng/mL, respectively.
- Badly hemolyzed serum or plasma samples should not be used, as they may give inaccurate results.
- Insufficient washing will result in poor precision and false absorbance readings.

**PERFORMANCE CHARACTERISTICS**

Parallelism:

Performed using a sequence of serial dilutions.



Standards

$y = 0.0505x + 0.0869$

Sample

$y = 0.0505x + 0.0869$

Recovery:

Sample was diluted 1:10 in 1x EIA Buffer prior to starting the assay. Three analyzed volumes of MPO were added to dispensed volumes of the diluted sample.

	% Recovery
Concentration 1	83
Concentration 2	90
Concentration 3	180

---

**REFERENCES**

---

1. Klebanoff, S. J.; (2005) *J of Leuk. Bio.* **77**:598
2. Klein, J.; (1982) *Immunology The Science of Self-Nonself Discrimination* p 427
3. Renliang, Z., *et al.*; (2001) *JAMA* **286(17)**:2136, 2140
4. Reumaux, D., *et al.*; (2003) *J of Leuk. Bio.* **73**:846
5. Van der Woode, F.J., and Schmitt, W.H.; (2004) Clinical Applications of Anti-neutrophil Cytoplasmic Antibody Testing. <http://www.medscape.com/viewarticle/467031>. Accessed July 31, 2008

---

**DISCLAIMER**

---

This information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Oxford Biomedical Research, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See catalog for additional terms and conditions of sale.

---

**ORDERING INFORMATION**

---

For additional kits or a complete catalog please call 800-692-4633.

---

**TECHNICAL SUPPORT**

---

If you need technical information or assistance with assay procedures, call our Technical Support Department at 800-692-4633 or 248-852-8815. Our staff will be happy to answer your questions about this or any other product in the Oxford Biomedical line.

---

**GUARANTEE AND LIMITATION OF REMEDY**

---

Oxford Biomedical Research, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the materials in this kit, except that these materials and this kit will meet our specifications at the time of delivery. Buyer's remedy and Oxford Biomedical Research, Inc.'s sole liability hereunder is limited to, at Oxford Biomedical Research, Inc.'s option, refund of the purchase price of, or the replacement of, material that does not meet our specification. By acceptance of our products, Buyer indemnifies and holds Oxford Biomedical Research, Inc. harmless against, assumes all liability for the consequence of its use or misuse by the Buyer, its employees, or others. Said refund or replacement is conditioned on Buyer notifying Oxford Biomedical Research, Inc. within thirty (30) days of the receipt of product. Failure of Buyer to give said notice within thirty (30) days of receipt of product shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

Oxford Biomedical Research, Inc.  
P.O. Box 522  
Oxford, MI 48371 U.S.A.

Orders: 800-692-4633  
Technical Service: 248-852-8815  
Fax: 248-852-4466  
E-mail: [info@oxfordbiomed.com](mailto:info@oxfordbiomed.com)

Made in the U.S.A.