

Total Antigen ELISA for Rat Prorenin/Renin

For Research Use Only

INTRODUCTION

Prorenin is a glycosylated aspartic protease that consists of 2 homologous lobes and is the precursor of renin. Renin activates the renin-angiotensin system by cleaving angiotensinogen, produced by the liver, to yield angiotensin I, which is further converted into angiotensin II by ACE, the angiotensin-converting enzyme primarily within the capillaries of the lungs. It has been reported that the levels of circulating prorenin (but not renin) are increased in diabetic subjects (1).

PRINCIPLES OF PROCEDURE

This ELISA uses a capture antibody coated to the 96-well plate to bind Rat Prorenin and Renin. A detection antibody conjugated to biotin is then applied. After washing, avidin conjugated to horseradish peroxidase (HRP) is applied. TMB substrate is added for color development which is proportional to the total concentration of prorenin and rennin in the sample. Sample concentrations can be determined by comparing OD values to the standard curve.

MATERIALS PROVIDED

Component	Contents	Quantity	Storage	Cat. No.
Coated Plate	Anti-Rat Prorenin/Renin	1 vial	4°C	RN46a
Standard	Rat Prorenin (lyophilized)	1 vial	4°C	RN46b
Wash Buffer	10x solution for washing plate	50 mL	4°C	RN46c
Primary Antibody	Anti Rat Prorenin/Renin Antibody	1 vial	4°C	RN46d
Avidin HRP	HRP Labeled Avidin	1 vial	4°C	RN46e
Substrate	TMB Substrate	10 mL	4°C	RN46f

MATERIALS NEEDED BUT NOT PROVIDED

1. Pipettes covering 0-10 μ l and 200-1000 μ l and tips
2. 12-channel pipette covering 30-300 μ l
3. 1N Sulfuric Acid (H₂SO₄)
4. Deionized Water (DI water)
5. Microtiter plate spectrophotometer with a 450 nm filter
6. Microtiter plate shaker with uniform horizontally circular movement up to 300 rpm

STORAGE CONDITIONS

1. Store this kit and its components at 4°C until use.
2. Reconstituted standards and primary may be stored at -70°C for later use. **DO NOT** freeze/thaw the Standards and Primary Antibody more than once.

PROCEDURAL NOTES

1. Use aseptic technique when opening and dispensing reagents.
2. This kit is designed to work properly as provided and instructed. Additions, deletions, or substitutions to the procedure or reagents are not recommended, as they may be detrimental to the assay.
3. Exercise universal precautions during the performance or handling of this kit or any component contained therein.

SAMPLE COLLECTION AND PREPARATION

Typical rat prorenin levels range from 0 to 400 ng/mL. The assay range for this ELISA is 0.1 to 100 ng/mL. Samples with prorenin levels above 100 ng/mL should be diluted in Blocking Buffer and retested.

REAGENT PREPARATION

The following solutions should be prepared fresh before starting the assay.

1. **TBS Buffer:** 0.1 M TRIS, 0.15 M NaCl, pH 7.4
2. **3% BSA Blocking Buffer:** 3% BSA in TBS Buffer.
3. **10x Wash Buffer:** Dilute the 50 mL of concentrate to 1x with 450 mL of DI water prior to use.
4. **Standard:** Reconstitute with 1.0 mL of 3% BSA Blocking Buffer and vortex gently to mix. Prepare according to included Standard Dilution Table immediately prior to use.
5. **Primary Antibody:** Reconstitute with 10 mL of 3% BSA Blocking Buffer and vortex gently to mix. Prepare immediately prior to use.
6. **Avidin HRP:** Dilute with 10 mL of 3% BSA Blocking Buffer and vortex gently to mix. Prepare immediately prior to use.

STANDARD PREPARATION

Reconstitute the Standard as directed on the vial to give a 1000 ng/mL Standard Stock Solution (see Reagent Preparation). **Do not prepare the standards until you are ready to apply them to the plate.**

Table 1: Preparation of Standard Curve

Standard	Prorenin Concentration (ng/mL)	Blocking Buffer (μ L)	Transfer Volume (μ L)	Transfer Source	Final Volume (μ L)
S ₁₀	100	900	100	Stock	750
S ₉	50	250	250	S ₁₀	300
S ₈	20	300	200	S ₉	250
S ₇	10	250	250	S ₈	250
S ₆	5	250	250	S ₇	300
S ₅	2	300	200	S ₆	250
S ₄	1	250	250	S ₅	250
S ₃	0.5	250	250	S ₄	300
S ₂	0.2	300	200	S ₃	250
S ₁	0.1	250	250	S ₂	500
S ₀	0	500	---	---	500

ASSAY PROCEDURE

1. Add 100 μ l of the Standards and unknowns to the wells in duplicate. Shake the plate at 300 rpm for 30 minutes at room temperature (RT). For a suggested plate layout, see **Scheme I** below.
2. Wash the plate 3 times according to the following wash procedure:
 - a. Remove the contents of each well by inversion of the plate.
 - b. Tap out the remaining contents of the plate onto a lint free paper towel.
 - c. Add 300 μ L of 1x Wash Buffer.
 - d. Let stand for 2-3 minutes.
 - e. Repeat procedure two more times, then proceed to step “F”.
 - f. Remove the contents of each well by inversion of plate into an appropriate disposal device.
 - g. Tap out the remaining contents of the plate onto a lint free paper towel, then proceed to step 3.
3. Add 100 μ l of the Primary Antibody to each well. Shake the plate at 300 rpm for 30 minutes at RT.
4. Wash the plate three times as in step 2.
5. Add 100 μ l of the Secondary Antibody to each well. Shake the plate at 300rpm for 30 minutes at RT.
6. Wash the plate three times as in step 2.
7. Add 100 μ l of TMB Substrate to each well. Shake the plate at 300 rpm for 2-10 minutes at RT.
8. Stop the reaction by adding 50 μ l of 1N H₂SO₄ to each well and read the plate at 450 nm.


Scheme I: Sample Plate Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	S10	S9	S8	S7	S6	S5	S4	S3	S2	S1	S0	U1
B	S10	S9	S8	S7	S6	S5	S4	S3	S2	S1	S0	U1
C	U2	U3	U4	U5	U6	U7	U8	U9	U10	U11	U12	U13
D	U2	U3	U4	U5	U6	U7	U8	U9	U10	U11	U12	U13
E	U14	U15	U16	U17	U18	U19	U20	U21	U22	U23	U24	U25
F	U14	U15	U16	U17	U18	U19	U20	U21	U22	U23	U24	U25
G	U26	U27	U28	U29	U30	U31	U32	U33	U34	U35	U36	U37
H	U26	U27	U28	U29	U30	U31	U32	U33	U34	U35	U36	U37

CALCULATIONS

1. Plot the A₄₅₀ against the concentration of uPA in the standards.
2. Fit a straight line through the points using a linear fit procedure.
3. Calculate the uPA concentrations in the unknowns using the equation generated by the standard curve.

Figure 1: Typical Standard Curve



EXPECTED VALUES

Rat prorenin levels range from 0-400 ng/ml depending on assay methodology². Human plasma levels of prorenin are greater in males than females and correlate positively with age and negatively with blood pressure³. Plasma and serum concentrations increase in several conditions such as pregnancy, progressive diabetes mellitus, diabetes mellitus with microvascular disease, and diabetic retinopathy^{4,5}.

REFERENCES

1. Luetscher, J.A., *et. al.*; (1985) N. Engl. J. Med. **312**:1412-1417
2. Ioannou, P., *et. al.*; (1991) Can. J. Physiol. Pharmac. **69**(9):1331-1340
3. Danser, A.H., *et. al.*; (1988) J. Hypertens. **16**:853-862
4. Yokota, H., *et. al.*; (2005) Br. J. Ophthalmol. **89**:871-873
5. Schmieder, R.E.; (2007) J. Hypertens. **25**:1323-1326

DISCLAIMER

This information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Oxford Biomedical Research, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See catalog for additional terms and conditions of sale.

TECHNICAL SUPPORT

If you need technical information or assistance with assay procedures, please call our Technical Support Department at 1-800-692-4633 or 1-248-852-8815. Our staff will answer your questions about this or any other product in the Oxford Biomedical line.

GUARANTEE AND LIMITATION OF REMEDY

Oxford Biomedical Research, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the material in this ELISA kit, except that these materials and this kit will meet our specifications at the time of delivery. Buyer's remedy and Oxford Biomedical Research, Inc.'s sole liability hereunder is limited to, at Oxford Biomedical Research, Inc.'s option, refund of the purchase price of, or the replacement of, material that does not meet our specification. By acceptance of our products, Buyer indemnifies and holds Oxford Biomedical Research, Inc. harmless against, assumes all liability for the consequence of its use or misuse by the Buyer, its employees, or others. Said refund or replacement is conditioned of Buyer notifying Oxford Biomedical Research, Inc. within (30) days of the receipt of product. Failure of Buyer to give said notice within said thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

Oxford Biomedical Research, Inc.
P.O. Box 522
Oxford, MI 48371 U.S.A.

Orders: 800-692-4633
Technical Service: 248-852-8815
Fax: 248-852-4466
E-mail: info@oxfordbiomed.com

Made in the U.S.A.