

**Product Specification Sheet**

Description: Anti-6-keto-PGF1alpha  
Product Number: PG40  
Aliquot: 1 mL  
Lot Number: pg40.030597DM  
Storage: -70° C

<b>DESCRIPTION:</b>	Sheep polyclonal antiserum specific for 6-keto-PGF1 alpha in PBS with 0.1% sodium azide.																						
<b>MOLECULAR WEIGHT:</b>	Target MW is 370.5																						
<b>STORAGE BUFFER:</b>	PBS with 0.1% sodium azide.																						
<b>STORAGE:</b>	Store at -70°C																						
<b>IMMUNOGEN:</b>	6-keto-PGF1 alpha																						
<b>USE:</b>	Thaw antiserum at room temperature. For radioimmunoassay procedures, make a 1:25 working dilution in PBS plus 0.1% gelatin.																						
<b>VOLUME:</b>	1.0 mL. Sufficient for at least 100 tubes in conventional liquid phase radioimmunoassay procedures.																						
<b>SPECIFICITY:</b>	<table><tr><td>6-keto-PGF1a</td><td>100%</td></tr><tr><td>PGF2a</td><td>2.1%</td></tr><tr><td>Thromboxane B2</td><td>0.07%</td></tr><tr><td>13,14-dihydro-15-keto PGF2a</td><td>0.03%</td></tr><tr><td>13,14-dihydro-6,15-diketo PGF2a</td><td>0.02%</td></tr><tr><td>PGD2</td><td>0.03%</td></tr><tr><td>PGE1</td><td>0.06%</td></tr><tr><td>PGE2</td><td>0.05%</td></tr><tr><td>6-keto-PGE1</td><td>8.6%</td></tr><tr><td>PGF1a</td><td>7.9%</td></tr><tr><td>Arachidonic Acid</td><td>0.02%</td></tr></table>	6-keto-PGF1a	100%	PGF2a	2.1%	Thromboxane B2	0.07%	13,14-dihydro-15-keto PGF2a	0.03%	13,14-dihydro-6,15-diketo PGF2a	0.02%	PGD2	0.03%	PGE1	0.06%	PGE2	0.05%	6-keto-PGE1	8.6%	PGF1a	7.9%	Arachidonic Acid	0.02%
6-keto-PGF1a	100%																						
PGF2a	2.1%																						
Thromboxane B2	0.07%																						
13,14-dihydro-15-keto PGF2a	0.03%																						
13,14-dihydro-6,15-diketo PGF2a	0.02%																						
PGD2	0.03%																						
PGE1	0.06%																						
PGE2	0.05%																						
6-keto-PGE1	8.6%																						
PGF1a	7.9%																						
Arachidonic Acid	0.02%																						
<b>SPECIES SPECIFICITY:</b>	Not species specific.																						
<b>REFERENCES:</b>	J.L. Fahey, et al., Amer. Soc. Microbiol., 685-691 (1986) W.S. Powell, Methods in Enzymology, 86; 467-476																						