

Product Specification Sheet

Description: Human Factor X
Product Number: CF35
Aliquot: 0.5 mg
Lot Number: cf35.4215
Storage: -20°C

DESCRIPTION:	Human Factor X purified from human plasma
MOLECULAR WEIGHT:	58.9 kDa
CONCENTRATION:	7.6 mg/ml
STORAGE BUFFER:	50% Glycerol
STORAGE:	Store at -20°C
SPECIFIC ACTIVITY:	122 U/mg in a Factor X clotting assay
PURIFICATION:	Human factor X is purified using a combination of conventional techniques and immunoaffinity chromatography.
PURITY:	>95% bt SDS-PAGE analysis