



OXFORD BIOMEDICAL RESEARCH

### Product Specification Sheet

Anti-Thrombopoietin (Mouse) Biotin  
Labeled  
Product Number: CT11  
Aliquot: 1.0 mg  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal sheep antibody is specific for mouse thrombopoietin and is provided as a biotin labeled IgG fraction.
<b>MOLECULAR WEIGHT:</b>	N/A
<b>CONCENTRATION:</b>	10 mg/mL
<b>SIZE:</b>	1.0 mg
<b>VOLUME:</b>	0.1 mL
<b>STORAGE:</b>	-70°C – <b>Avoid repeated freeze/thaw cycles.</b>
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Sheep serum
<b>BUFFER:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>PURIFICATION:</b>	IgG fraction purified by immobilized Protein G.
<b>SPECIES SPECIFICITY:</b>	Mouse

For research purposes only.