



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Sheep Anti-Rabbit Fibrinogen IgG Fraction
Product Number: CT32
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Anti-Fibrinogen (Rabbit), Sheep Polyclonal Antibody, IgG Fraction
MOLECULAR WEIGHT:	160,000
CONCENTRATION:	10 mg/mL
SIZE:	1.0 mg
VOLUME:	0.1 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Sheep serum
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
PURIFICATION:	IgG fraction purified by immobilized Protein G

For research purposes only.