

Product Specification Sheet

Description: Anti-Myeloperoxidase, Human Neutrophil
Product Number: GP50
Aliquot: 1.0 mL
Lot Number: gp50.2559
Storage: -70°C

DESCRIPTION:	This is a polyclonal antibody, IgG fraction, was prepared in rabbit using purified, full length human neutrophil myeloperoxidase.
MOLECULAR WEIGHT:	59 kDa
CONCENTRATION:	1.01 mg/ml
STORAGE BUFFER:	PBS; 0.1% Azide
STORAGE:	Store at -70°C
IMMUNOGEN:	Human full-length native neutrophil myeloperoxidase
PURIFICATION:	Purified using immobilized Protein A
CLASS:	IgG
APPLICATIONS:	ELISA, WB
SPECIES SPECIFICITY:	Human. Has not been evaluated in other species.