



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Endothelial Nitric Oxide Synthase (eNOS)

Product Number: NS 15

Aliquot: 100  $\mu$ L

Lot Number: ns15.P9013A

Storage: -80°C

<b>DESCRIPTION:</b>	Anti-eNOS, Endothelial Nitric Oxide Synthase, rabbit polyclonal antibody. This antibody does not cross react with iNOS or bNOS.
<b>MOLECULAR WEIGHT:</b>	135 kDa
<b>CONCENTRATION:</b>	Not Determined
<b>PURIFICATION:</b>	None
<b>STORAGE BUFFER:</b>	PBS; 0.05% Sodium Azide
<b>STORAGE:</b>	-80°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	A specific peptide corresponding to acids 599-613 from bovine eNOS.
<b>APPLICATIONS:</b>	Western: 1:1000 Immunocytochemistry: 1:250
<b>SPECIES SPECIFICITY:</b>	This antibody is known to be reactive in human, mouse, rat, and bovine.