



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Protein Kinase C β II
Product Number: PK 20
Aliquot: 1.0 mL
Lot Number: 044693JC
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C β II (PKC β II).
MOLECULAR WEIGHT:	77 kDa Calculated
CONCENTRATION:	Not Determined
ISOTYPE:	NA
STORAGE BUFFER:	Phosphate Buffered Saline; 0.05% Sodium Azide
STORAGE:	-70°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
IMMUNOGEN:	KLH conjugated synthetic peptide corresponding to amino acids 661- 673 of PKC β II. Sequence: ⁶⁶¹ F-V-N-S-E-F-L-K-P-E-V-K-S ⁶⁷³
PURIFICATION:	NA
APPLICATIONS:	Western Blotting: 1:100 Dilution
SPECIES SPECIFICITY:	Human, mouse and rat among other mammalian species.
REFERENCE:	JBC (1993) 268 :20110-20115