



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA 800-692-4633 • Fax 248-852-4466

www.oxfordbiomed.com

Product specifications

Protein Kinase C Zeta (ζ) Synthetic Peptide

Product # PK 42 Typical Lot

Peptide Type:	Distinct region of protein kinase C C-terminalSequence: NH ₂ -C-F-E-G-F-E-Y-I-N-P-L-L-L-S-T-E-E-S-V-COOH
Applications:	Will compete with the anti-peptide antibody in Western blotting for confirmation of the PKC specificity of the antibody binding.
Contents per vial:	100 µg of peptide in 500 µL of ammonium bicarbonate buffer
Storage:	Stable at 4°C. For long term storage, aliquot product and freeze at -20°C. Avoid repeated freeze/thaw cycles.