



## OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA  
USA: 800-692-4633 • Fax: 248-852-4466  
[www.oxfordbiomed.com](http://www.oxfordbiomed.com)

### *Product specifications*

## **Anti-PKC Antibody Panel**

### **Product # PK 50**

The panel is useful for screening relative expression of PKC isoforms. Each vial contains sufficient antibody for at least 5 Western blots using Chemiluminescence.

This panel contains one vial each, of the following:

PK 10 – Anti-PKC alpha peptide antibody

PK 15 – Anti-PKC beta I peptide antibody

PK 20 – Anti-PKC beta II peptide antibody

PK 25 – Anti-PKC gamma peptide antibody

PK 30 – Anti-PKC delta peptide antibody

PK 35 – Anti-PKC epsilon peptide antibody

PK 40 – Anti-PKC zeta peptide antibody

PK 45 – Anti-PKC eta peptide antibody