



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466

www.oxfordbiomed.com

Product specifications

Polyclonal Protein Kinase Consensus Sequence Antibody

Product # PK 75 Typical Lot

- Immunoglobulin type:** Polyclonal rabbit IgG
- Specificity:** This antibody was generated against the G-T-P-E-Y-L-A-P-E peptide and is excellent for Western blot detection of protein kinases containing this sequence.
- Applications:** A dilution of 1: 50 is recommended for Western blotting using alkaline phosphatase. A further dilution may be possible with more sensitive detection systems.
- Note:** Use of 1% normalized goat serum added to a usual blocking agent is recommended to reduce non-specific interactions.
- Contents/vial:** 1.0 mL of rabbit antiserum in PBS
- Storage:** Stable for several weeks at 4°C, several months at -20°C, and indefinitely at -70°C. Do not refreeze this product.