

Product Specification Sheet

Description: Rabbit anti serine/threonine PP 2B α
Product Number: PP16
Aliquot: 100 uL
Lot Number: pp16.12095
Storage: -20C

DESCRIPTION:	Rabbit anti-serine/threonine protein phosphatase 2B α (anti-calcineurin)
TARGET MW:	60 kDa
CONCENTRATION:	1.0 mg/mL
STORAGE BUFFER:	PBS with 0.08% Sodium Azide
STORAGE:	Store at -20°C. Avoid freeze/thaw cycles.
IMMUNOGEN:	Purified Peptide NH ₂ -Arg-Met-Pro-Pro-Arg-Arg-Asp-Ala-Met-Pro-Ser-Asp-Ala-COOH (a.a. 482-494) conjugated to KLH.
PURIFICATION:	Ammonium Sulfate Precipitation
REACTIVITY:	Human, Mouse, Rat, Bovine
APPLICATIONS:	Western blotting and IHC