



## OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA  
USA: 800-692-4633 • Fax: 248-852-4466  
[www.oxfordbiomed.com](http://www.oxfordbiomed.com)

### *Product specifications*

## Anti Epoxide Hydrolase

### Product # PR 35      Typical Lot

**Immunoglobulin**

**Type:** Polyclonal goat IgG

**Specificity:**

Cross-reacts with rabbit and rat epoxide hydrolase (MW = 50 kDa). Non-specific cross-reactivity with 60 and 70 kDa species.

**Applications:**

Western blots

**Contents:**

100  $\mu$ L of goat IgG containing 0.05% sodium azide.

**Working Dilution:**

A dilution of 1/1000 for Western blots will yield a positive result with 0.2  $\mu$ g of purified rabbit or rat epoxide hydrolase using alkaline phosphatase conjugated donkey anti-goat antibodies according to the manufacturer's instructions.

**Storage:**

Store at -20°C. Avoid repeated freeze-thaw cycles.