

Mouse Anti-Human GST Z1-1
Product Number: GS21
Aliquot: 100 µg
Lot Number: gs21.150217
Storage: -80°C

DESCRIPTION:	This mouse monoclonal antibody to human glutathione S-transferase Z1-1 is protein A & A-G affinity purified and shows no cross reactivity for all 18 cytosolic GSTs.
CONCENTRATION:	1.00 mg/mL
STORAGE BUFFER:	Phosphate Buffered Saline, pH 7.4 and 25% Glycerol and 0.02% NaN ₃
STORAGE:	-80°C. Avoid repeated freeze thaw cycles
IMMUNOGEN:	Recombinant human GST Z1-1 expressed in <i>E. coli</i>
PURIFICATION:	Protein A & A-G affinity purified
APPLICATIONS:	Enzyme Immunoassay. Immunohistochemistry. Western Blotting: 1:4K dilution is reactive with ≥ 0.2µg of recombinant human GST Z1-1 (OBR product #GS77) in conjunction with an HRP conjugated anti-mouse secondary antibody.
SPECIES SPECIFICITY:	Human

