

Mouse Anti-Human GST O1-1  
Product Number: GS26  
Aliquot: 100 µg  
Lot Number: gs26.150205  
Storage: -80°C

<b>DESCRIPTION:</b>	This mouse monoclonal antibody to human glutathione S-transferase O1-1 is protein A & A-G affinity purified and shows no cross reactivity for all 18 cytosolic GSTs.
<b>CONCENTRATION:</b>	1.00 mg/mL
<b>STORAGE BUFFER:</b>	Phosphate Buffered Saline, pH 7.4 and 25% Glycerol and 0.02% NaN <sub>3</sub>
<b>STORAGE:</b>	-80°C. Avoid repeated freeze thaw cycles
<b>IMMUNOGEN:</b>	Recombinant human GST O1-1 expressed in <i>E. coli</i>
<b>PURIFICATION:</b>	Protein A & A-G affinity purified
<b>APPLICATIONS:</b>	Enzyme Immunoassay. Immunohistochemistry. Western Blotting: 1:5K dilution is reactive with ≥ 0.2µg of recombinant human GST O1-1 (OBR product #GS75) in conjunction with an HRP conjugated anti-mouse secondary antibody.
<b>SPECIES SPECIFICITY:</b>	Human

