

Product Specification Sheet

Calf Thymus Core Histones

Product Number: NF20

Aliquot: 100 µg

Lot Number: nf20.220302

Storage: -20°C

DESCRIPTION: 100 µg lyophilized, reconstitute in 1mL of DI Water or Buffer for a 100µg/mL solution.

STORAGE: -20°C; **AVOID MULTIPLE FREEZE-THAW CYCLES.**

PURITY: >90%

MOLECULAR WEIGHT: Mixture

SOURCE: Calf Thymus