

Product Specification Sheet

Anti-Factor IX (Human)
Antiserum
Product Number: CF02
Aliquot: 1.0 mL
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	This polyclonal antibody to human factor IX is produced in sheep and provided as antiserum.
MOLECULAR WEIGHT:	N/A
CONCENTRATION:	Whole Antiserum
VOLUME:	1.0 mL
STORAGE:	-70°C – Avoid repeated freeze/thaw cycles.
FORM:	Frozen Liquid
SOURCE:	Sheep Serum
SPECIES SPECIFICITY:	Human

For research purposes only.