

Product Specification Sheet

Dog Fibrinogen
Product Number: CT29
Aliquot: 1.0 mg
Lot Number: ct29.313
Storage: -70°C

DESCRIPTION:	Prepared from fresh dog plasma using several chromatographic steps.
MOLECULAR WEIGHT:	330,000
CONCENTRATION:	20 mg/mL
SIZE:	1.0 mg
VOLUME:	0.05 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SPECIFIC ACTIVITY:	96% clottable
PURITY:	>95% by SDS-PAGE analysis
SOURCE:	Dog plasma
BUFFER:	0.02M Sodium Citrate-HCl; pH 7.4
THAWING:	Thaw at 37°C without agitation until completely liquid. Mix gently by pipet prior to use. Agitation may cause clotting. Keep at 25-37°C during use. Make aliquots and flash freeze unused portion.

For Research Purposes Only.