



Product Specification Sheet

Description: Human Alpha-Thrombin
Product Number: CT36
Aliquot: 1.0 mg
Lot Number: ct36.bb0926
Storage: -20°C

DESCRIPTION:	Human alpha thrombin prepared from purified prothrombin obtained from human plasma.
MOLECULAR WEIGHT:	36.7 kDa
CONCENTRATION:	7.7 mg/mL
STORAGE BUFFER:	50% glycerol
STORAGE:	-20°C
SPECIFIC ACTIVITY:	3367 U/mg as determined by a fibrinogen clotting assay
PURITY:	>95% by SDS-PAGE analysis