



## Product Specification Sheet

Description: Rabbit Anti-human Fibrinogen, IgG Fraction  
Product Number: CT45  
Aliquot: 1.0 mg  
Lot Number: ct45.609  
Storage: -70°C

<b>DESCRIPTION:</b>	Rabbit Anti-human Fibrinogen, IgG Fraction
<b>MOLECULAR WEIGHT:</b>	340,000
<b>CONCENTRATION:</b>	10.0 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Full length native human fibrinogen
<b>PURIFICATION:</b>	IgG fraction by immobilized Protein A
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	ELISA, WB
<b>SPECIES SPECIFICITY:</b>	Human