

Product Specification Sheet

Goat Anti-Dinitrophenol (DNP)
IgG
Product Number: D04
Aliquot: 1.0 mg
Lot Number: d04.
Storage: -70°C

DESCRIPTION:	This goat polyclonal antibody is specific to DNP and is intended for use in experiments related to receptor-mediated endocytosis and oxidative stress.
CONCENTRATION:	10 mg/mL as determined by BCA assay using immunoglobulin as a standard.
STORAGE BUFFER:	PBS
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	DNP conjugated to KLH
PURIFICATION:	Ammonium Sulfate Fractionation
APPLICATIONS:	This antibody has been shown to have utility in ELISA and Western Blot. Other applications have not been evaluated.
SPECIES SPECIFICITY:	Not species specific. May not exhibit utility when using in a goat experimental model.