

Product Specification Sheet

Human Multimeric Vitronectin
Product Number: GP02
Aliquot: 100 µg
Lot Number: gp02.605A
Storage: -70°C

DESCRIPTION:	This human multimeric vitronectin is urea prepared from human plasma using non-denaturing chromatography.
MOLECULAR WEIGHT:	75 kDa
CONCENTRATION:	3.4 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; pH 7.4
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.