

**Product Specification Sheet**

Mouse Vitronectin  
Product Number: GP04  
Aliquot: 0.1 mg  
Lot Number: gp04.808  
Storage: -70°C

|                          |  |
|--------------------------|--|
| <b>DESCRIPTION:</b>      | Mouse Vitronectin, Urea Prepared             |
| <b>ALIQUOT:</b>          | 0.1 mg in 0.042 mL                           |
| <b>CONCENTRATION:</b>    | 2.4 mg/mL                                    |
| <b>MOLECULAR WEIGHT:</b> | ~75 KDa                                      |
| <b>STORAGE BUFFER:</b>   | 0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.4 |
| <b>STORAGE:</b>          | Store frozen at -70°C                        |