

Product Specification Sheet

Rat Vitronectin (Multimeric)
Product Number: GP05
Aliquot: 100 µg
Lot Number: gp05.313
Storage: -70°C

DESCRIPTION: This rat multimeric vitronectin is urea prepared from rat plasma.

MOLECULAR WEIGHT: 75 kDa

CONCENTRATION: 0.95 mg/mL determined spectrophotometrically.

STORAGE BUFFER: 0.05 M Sodium Phosphate; 0.1 M NaCl; pH 7.4

STORAGE: -70°C. **AVOID REPEATED FREEZE THAWS.**

PURITY: > 95% pure by SDS-PAGE