

Product Specification Sheet

Rabbit Anti-Mouse Vitronectin
(IgG Fraction)
Product Number: GP16
Aliquot: 1.0 mg
Lot Number: gp16.105
Storage: -70°C

DESCRIPTION:	Rabbit polyclonal antibody against vitronectin.
MOLECULAR WEIGHT:	75 kDa
CONCENTRATION:	3.4 mg/mL
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Vitronectin
PURIFICATION:	The IgG fraction was purified by immobilized Protein A.
APPLICATIONS:	ELISA, IHC, Western Blot
SPECIES SPECIFICITY:	Specific for mouse.