

Product Specification Sheet

Description: Anti-Rat GST Yb (Mu)
Product Number: GS23
Aliquot: 100uL
Lot Number: gs23.151116
Storage: -70°C

DESCRIPTION:	Goat antiserum against Rat GST Yb (Mu class)
CONCENTRATION:	68 mg/mL Neat Antiserum
STORAGE:	Store at -70°C
IMMUNOGEN:	Affinity purified Rat GST Yb1
REFERENCES:	Holmes, W., et al., Br. J. Cancer, 62: 209-212 (1990)
DILUTION:	A dilution of 1:1000 will yield a positive result with 0.1 ug of purified rat Yb1 GST or recombinant human M1-1.
APPLICATIONS:	Western blotting
SPECIES SPECIFICITY:	Rat and Human GST Mu