

Product Specification Sheet

Anti-15-Lipoxygenase Form-2
Product Number: LX25
Aliquot: 100 µL
Lot Number: 180103DM
Storage: -70°C

DESCRIPTION:	15-Lipoxygenase Form-2 (15-LOX-2) is expressed primarily in lung tissue, hair roots, the prostate and the cornea and is a useful biomarker in eicosanoid studies.
MOLECULAR WEIGHT:	78.8 kDa
CONCENTRATION:	Not Determined
PRESENTATION:	Rabbit Antisera
STORAGE:	-70°C
APPLICATIONS:	This antibody has been tested in Western blot, immunohistochemistry (paraffin tissues), immunoprecipitation and ELISA applications. Recommended Western dilution: 1:1,000
IMMUNOGEN:	Recombinant Human 15-Lipoxygenase-form-2.
SPECIES SPECIFICITY:	This antibody will detect 15-LOX-2 in Human and 8S-LOX in mouse.
REFERENCE:	Jisaka, M. et. al. (1997) J. of Biol. Chem. 272: 39 pp.24410-24416.