

Product Specification Sheet

Mouse Anti-Human GST Z1-1
 Product Number: GS21
 Aliquot: 100 µg
 Lot Number: gs21.150217
 Storage: -70°C

- DESCRIPTION:** This mouse monoclonal antibody to human glutathione S-transferase Z1-1 is protein A & A-G affinity purified and shows no cross reactivity for all 18 cytosolic GSTs.
- CONCENTRATION:** 1.00 mg/mL
- STORAGE BUFFER:** Phosphate Buffered Saline, pH 7.4 and 25% Glycerol and 0.02% NaN₃
- STORAGE:** -70°C. Avoid repeated freeze thaw cycles
- IMMUNOGEN:** Recombinant human GST Z1-1 expressed in *E. coli*
- PURIFICATION:** Protein A & A-G affinity purified
- APPLICATIONS:** Enzyme Immunoassay. Immunohistochemistry. Western Blotting: 1:4K dilution is reactive with ≥ 0.2µg of recombinant human GST Z1-1 (OBR product #GS77) in conjunction with an HRP conjugated anti-mouse secondary antibody.
- SPECIES SPECIFICITY:** Human

