

**Product Specification Sheet**

Human tPA, Immobilized  
Product Number: PA48  
Aliquot: 1.0 mg  
Lot Number: Typical Lot  
Storage: 4°C

<b>DESCRIPTION:</b>	May be used to immunopurify monoclonal and polyclonal antibodies directed against human tPA. May be used repeatedly
<b>MOLECULAR WEIGHT:</b>	70,000
<b>CONCENTRATION:</b>	2.0 mg/mL coupled
<b>SIZE:</b>	1.0 mg
<b>VOLUME:</b>	2.0 mL of a 1:1 slurry
<b>STORAGE:</b>	4°C
<b>FORM:</b>	Resin
<b>SOURCE:</b>	CHO cells
<b>BUFFER:</b>	0.1M HEPES; 0.1M NaCl; 0.02% Sodium Azide; pH 7.4

For Research Purposes Only.