



## Product Specification Sheet

Recombinant Mouse tPA -  
Active  
Product Number: PA56  
Aliquot: 100  $\mu$ g  
Lot Number pa56.618  
Storage: -70°C

<b>DESCRIPTION:</b>	This mouse tPA is enzymatically active and is expressed from baculoviral containing insect cells.
<b>MOLECULAR WEIGHT:</b>	70 kDa
<b>CONCENTRATION:</b>	1.0 mg/mL
<b>STORAGE BUFFER:</b>	0.4 M HEPES; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C
<b>ACTIVITY:</b>	>95% active
<b>PURITY:</b>	>95% pure as determined by SDS-PAGE analysis
<b>SOURCE:</b>	Insect cell culture