


Anti-Human tPA MoAB
Product Number: PA64
Aliquot: 1.0 mg
Lot Number: pa64.323
Storage: -70°C

DESCRIPTION:	This is a monoclonal antibody produced in mouse and reacts with both free and complex forms of tPA.
MOLECULAR WEIGHT:	70 kDa
CLONE:	629-2A153
CONCENTRATION:	1.96 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1mM NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant full length wild-type Human tPA (glycosylated)
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
PURITY:	≥ 98% assessed by visual inspection of a Coomassie [®] blue-stained SDS-PAGE gel.
APPLICATIONS:	Suitable for use in ELISA and Western blot.
SPECIES SPECIFICITY:	Human
CROSS-REACTIVITY:	Pig



1 2 3 4

629-2A153

1. Human tPA.
2. Human tPA in complex with PAI-1.
3. Retavase[®]Boehringer/Mannheim (tPA fragment).
4. Retavase[®] in complex with PAI-1.