



### Product Specification Sheet

Description: Preimmune Control Serum for Polyclonal Anti-PGHS-1  
Product Number: PG20c  
Aliquot: 1 mL  
Lot Number: 20c.041201SP  
Storage: -20°C

**APPLICATIONS:** Use at same dilution as PG20 for negative control with immunoassay or western blots.

**CONTENTS PER VIAL:** 1 mL of preimmune rabbit anti-serum.

**STABILITY:** Stable for several days at 4°C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freeze/thaw cycles.