

Product Specification Sheet

Description: Ac-TVASSSTA
Octapeptide
Product Number: PI14
Aliquot: 5 mg
Lot Number: pi14.m1660
Storage: -70°C

DESCRIPTION: Ac-TVASSSTA Octapeptide mimics the N-terminal residues of the reactive center loop of PAI-1. This octapeptide will form a stable complex converting PAI-1 into a substrate for tPA.

MOLECULAR WEIGHT: 765

AMOUNT: 5 mg of lyophilized powder

STORAGE: -70°C

PURITY: >95% by HPLC analysis