

Product Specification Sheet

Description: Murine PAI-1
(Stable Mutant)
Product Number: PI18
Aliquot: 0.5 mg
Lot Number: pi18.203
Storage: -70°C

DESCRIPTION: Murine Plasminogen Activator Inhibitor-1 (Stable Mutant)
produced in *E. coli*

MOLECULAR WEIGHT: 43 kDa

CONCENTRATION: 2.3 mg/mL

STORAGE BUFFER: 0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0

STORAGE: -70°C

PURITY: >95% by SDS-PAGE analysis