

### Product Specification Sheet

Description: Anti-Protein Kinase C Beta II  
Product Number: PK23  
Aliquot: 100 uL  
Lot Number: pk23.101195yy  
Storage: -20°C

<b>DESCRIPTION:</b>	Anti-Protein Kinase C Beta II polyclonal antibody
<b>MOLECULAR WEIGHT:</b>	77 kDa
<b>CONCENTRATION:</b>	Antiserum
<b>STORAGE BUFFER:</b>	Antiserum containing 0.05% sodium azide
<b>STORAGE:</b>	Stable for several months at 4°C, at least 1 year at -20°C and indefinitely at -70°. Avoid repeated freeze/thaw cycles.
<b>IMMUNOGEN:</b>	Prepared against synthetic peptide corresponding to a distinct C terminal region, (SFVNSEFLKPEVKS, residues 660-673) from rat protein kinase C $\beta$ II. Based on known homologies it is expected to react with most mammalian species.
<b>APPLICATIONS:</b>	Western blotting, ELISA, IHC
<b>REFERENCES:</b>	1. Wetsel, WC, et al., J. Cell Biol., 117:21-33 (1992)