

Protein Kinase C Delta (δ) Standard

Product # : PK 31

Lot # : 001993JC

Contents/vial:	1.0 μ g of protein in 50 μ L.
Description:	Human recombinant kinase C δ expressed in baculoviral insect expression system. The protein is suspended in SDS/mercaptoethanol sample buffer with a bromphenol blue dye marker.
Protein Concentration:	20 μ g/mL
Enzyme activity:	Not enzymatically active
Applications:	Intended for use as a standard for Western blots using enhanced chemiluminescence. For blots using less sensitive substrates, such as alkaline phosphatase, larger amounts of protein should be loaded onto the gel.
Recommended Dilution:	Load 2 - 5 μ L of standard directly from vial. No dilution is necessary.
Storage:	Stable indefinitely at -70°C; for at least one year at -20°C.