

Product Specification Sheet

Description: Anti-Protein Kinase C Zeta
Product Number: PK43
Aliquot: 100 uL
Lot Number: pk43.033195SZ
Storage: -20°C

DESCRIPTION:	Anti-Protein Kinase C Zeta polyclonal antibody
MOLECULAR WEIGHT:	71 kDa
CONCENTRATION:	Antiserum
STORAGE BUFFER:	Antiserum containing 0.05% sodium azide
STORAGE:	Stable for several months at 4°C, up to a year at -20°C, and indefinitely at -70° C. Avoid repeated freeze thaws.
IMMUNOGEN:	Prepared against synthetic peptide (⁵⁷⁹ EYINPLLLSAEESV ⁵⁹²) corresponding to a distinct C terminal region from rat protein kinase C ζ and based on known homologies is expected to react with the PKC ortholog in most mammalian species.
APPLICATIONS:	Western blotting, IHC, ELISA
REFERENCES:	1. Wetsel, W. C., et al., <i>J. Cell Biol.</i> , 117 : 21-133 (1992)