

Product Specification Sheet

Description: Anti-PKC Antibody Panel
Product Number: PK50
Aliquot: 1 vial each
Lot Number: Per individual vial
Storage: -20°C

DESCRIPTION:

This panel is useful for screening relative expression of PKC isoforms. Each vial contains sufficient antibody for at least 5 Western blots using chemiluminescence detection.

This panel contains one vial each, of the following:

- PK10 – Anti-PKC alpha peptide antibody
- PK15 – Anti-PKC beta I peptide antibody
- PK20 – Anti-PKC beta II peptide antibody
- PK25 – Anti-PKC gamma peptide antibody
- PK30 – Anti-PKC delta peptide antibody
- PK35 – Anti-PKC epsilon peptide antibody
- PK40 – Anti-PKC zeta peptide antibody
- PK45 – Anti-PKC eta peptide antibody