

Product Specification Sheet

Description: Mouse Plasminogen
Product Number: PL28
Aliquot: 1.0 mg
Lot Number: pl28.311
Storage: -70°C

DESCRIPTION:	Mouse plasminogen prepared from fresh mouse plasma
MOLECULAR WEIGHT:	92 kDa
CONCENTRATION:	5.9 mg/ml
STORAGE BUFFER:	0.1 M HEPES; 0.1 M NaCl; pH 7.4
STORAGE:	Store up to 3 years at -70°C; avoid multiple freeze/thaw cycles
PURIFICATION:	Immobilized lysine chromatography
PURITY:	>95% by SDS-PAGE analysis