

Product Specification Sheet

Sheep anti-Human Plasminogen
Polyclonal Antibody
Product Number: PL70
Aliquot: 1.0 mg/mL
Lot Number: pl70.507
Storage: -70°C

DESCRIPTION:	Polyclonal Sheep anti-Human Plasminogen IgG Fraction
MOLECULAR WEIGHT:	92,000
CONCENTRATION:	5.1 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate, 0.1M NaCl, 1mM EDTA, pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Purified Human full length native protein.
PURIFICATION:	IgG fraction purified by immobilized Protein G.
CLASS:	IgG
APPLICATIONS:	ELISA, Western Blotting
SPECIES SPECIFICITY:	Human