

Product Specification Sheet

Anti-Rat CYP450 2E1 MoAB
Product Number: PM32
Aliquot: 100 µL
Lot Number 040700SV
Storage: -70°C

DESCRIPTION:	This mouse monoclonal antibody is specific for rat cytochrome P450 2E1 and is provided as purified IgG from ascites fluid.
CONCENTRATION:	Not Determined
STORAGE:	-70°C. AVOID MULTIPLE FREEZE THAW CYCLES
IMMUNOGEN:	Recombinant rat cytochrome P450 2E1
APPLICATIONS:	Western: 1:1000
MOLECULAR WEIGHT:	58.6 kDa, Calculated
CLONE NUMBER:	1-98-1
ISOTYPE:	Presented as IgG Purified MoAB from ascites.
NOTES:	This AB has been tested with 2.5 µg and 10 µg of pyridine induced rat liver microsomes using ECL and alkaline phosphatase detection methods respectively. Membrane overexposure to this AB may result in cross reactivity with other P450 isoforms.
REFERENCES:	<ol style="list-style-type: none">1. Gelboin et al. (1985) <i>J. Biochem. Pharmacol.</i> 34:2225-2234.2. Gelboin et al. (1993) <i>J. Biochem. Pharmacol.</i> 46:787-790.