

### Product Specification Sheet

Anti-Rabbit Cytochrome b5  
Product Number: PR91  
Aliquot: 100 µL  
Lot Number pr91.070712  
Storage: -70°C

<b>DESCRIPTION:</b>	This is a goat polyclonal antibody and is specific for cytochrome b5 in rabbit but cross reacts with human and rat. This antibody is non-cross reactive with the 50 and 60 kDa species of cytochrome b5.
<b>CONCENTRATION:</b>	Not determined.
<b>PRESENTATION:</b>	Neat antiserum with 0.08% (w/v) sodium azide.
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Rabbit cytochrome b5
<b>PURIFICATION:</b>	Not purified
<b>APPLICATIONS:</b>	This antibody is suitable for use in Western blotting at 1:1000 dilution.
<b>SPECIES SPECIFICITY:</b>	Rabbit, Rat and Human.