

Product Specification Sheet

Receptor Associated Protein –
Low Endotoxin
Product Number: RA02
Aliquot: 500 µg
Lot Number ra02.1019
Storage: -70°C

DESCRIPTION:	This human receptor associated protein (RAP) was made as a GST fusion protein and cleaved with thrombin. This RAP binds with high affinity to certain low-density lipoprotein (LDL) receptors such as LRP, gp330/megalin and the very low-density lipoprotein (VLDL).
MOLECULAR WEIGHT:	39 kDa
CONCENTRATION:	0.43 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.0067 M Phosphate; 0.15 M NaCl; pH 7.40
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
PURITY:	>95% pure by visual inspection of Coomassie™ Blue stained SDS-PAGE gel.
SOURCE:	E. coli