

Product Specification Sheet

Anti-Human Receptor Associate
Protein MoAB
Product Number: RA20
Aliquot: 1.0 mg
Lot Number ra20.807
Storage: -70°C

DESCRIPTION:	This is a mouse monoclonal antibody is specific for the 38 kDa low-density lipoprotein receptor-associated protein (RAP).
MOLECULAR WEIGHT:	38 kDa
CLONE:	MA7F1
CONCENTRATION:	5.3 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Tris-HCl; 0.1 M NaCl; pH 7.4
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS
IMMUNOGEN:	Human RAP
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
PURITY:	≥ 98% assessed by visual inspection of a Coomassie® blue-stained SDS-PAGE gel.
APPLICATIONS:	Suitable for use in Western blot at 1 µg/mL.
SPECIES SPECIFICITY:	Human. Does not cross-react with mouse or rabbit.