

Description: Anti-Human VLDL Receptor Monoclonal Antibody
Product Number: RA34
Aliquot: 1.0 mg
Lot Number: ra34.1008
Storage: -70°C

DESCRIPTION:	Anti-Human VLDL Receptor Monoclonal Antibody. Binds to an epitope on the ligand binding domain of the VLDL Receptor (AA 31-355).
MOLECULAR WEIGHT:	85 kDa
CONCENTRATION:	4.6 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate, 0.1M NaCl, 1mM EDTA, pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant fragment of soluble VLDL receptor containing ligand binding repeats 1-8.
PURIFICATION:	Immobilized protein A
CLASS:	IgG1
APPLICATIONS:	ELISA, Western Blotting
SPECIES SPECIFICITY:	Human. Cross reactivity to other species has not been tested.