

Product Specification Sheet

Description: Inhibitory Anti-human prorenin monoclonal antibody

Product Number: RN39

Aliquot: 0.1 mg

Lot Number: m39.913

Storage: -70°C

DESCRIPTION:	This mouse monoclonal antibody against human prorenin is specific to prorenin and does not bind or interact with active renin. The antibody is inhibitory and binds in the AA region of 30-43 of prorenin. It does not react with rat or mouse prorenin.
TARGET MW:	43.7 kDa
CONCENTRATION:	8.2 mg/ml
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	Store at -70°C
IMMUNOGEN:	Recombinant full length, wild type, glycosylated human prorenin
PURIFICATION:	IgG fraction
CLASS:	IgG2b kappa
APPLICATIONS:	WB, ELISA
SPECIES SPECIFICITY:	Human