

### Product Specification Sheet

Description: Human PAI-1 tPA complex  
Product Number: UP44  
Aliquot: 0.1 mg  
Lot Number: up44.120  
Storage: -70°C

<b>DESCRIPTION:</b>	This PAI-1 tPA complex was created by reacting a human PAI-1 stable mutant with an excess of human sc-tPA. Excess tPA was removed by immobilized PAI-1.
<b>MOLECULAR WEIGHT:</b>	113 kDa
<b>CONCENTRATION:</b>	1.0 mg/mL
<b>STORAGE BUFFER:</b>	0.4 M HEPES; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	Store at -70°C.
<b>SOURCE:</b>	PAI-1 mutant was produced in e. coli. tPA was produced in CHO cells.
<b>PURITY:</b>	>95% by SDS PAGE analysis.