

Product Specification Sheet

Description: Recombinant Human
HMW Urokinase
Product Number: VA50
Aliquot: 0.1 mg
Lot Number: va50.1016
Storage: -70°C

DESCRIPTION:	Recombinant Human HMW Urokinase produced in insect cell culture as the two chain form.
MOLECULAR WEIGHT:	54 kDa
CONCENTRATION:	1.7 mg/ml
STORAGE BUFFER:	0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0
STORAGE:	Store at -70°C
PURITY:	>95% by SDS-PAGE analysis
APPLICATIONS:	Plasminogen activation and receptor binding studies.